


# Surgilast®

## Tubular Elastic Dressing Retainer Filet tubulaire élastique à pansements

**23 meters** (25 yds. stretched)

**23 mètres** (25 verges à l'étirement)



**Glenwood Laboratories**  
Oakville, Ontario L6L 5M2

Surgilast® is a registered trademark of Glenwood.  
Surgilast® est une marque déposée de Glenwood.

Surgilast, a tubular elastic bandage retainer, is a stretchable net material used to hold dressings in place securely on all parts of the body without the use of adhesive tape. It is made from high quality nylon and rubber. Available in fifteen sizes, it is packaged in a dispenser box ready for use.

Le filet tubulaire élastique de maintien de pansements Surgilast est un matériel étirable spécialement conçu pour retenir les pansements bien en place, sans diachylon, quelle que soit la partie du corps. Il est fabriqué en nylon et caoutchouc de la plus grande qualité. Disponible en 15 tailles prêtes à l'emploi il est présenté en boîte distributrice.

| Cat. # | Size | Application                           |
|--------|------|---------------------------------------|
| GL-701 | 1    | Fingers, Toes, Wrist                  |
| GL-702 | 2    | Small Hand, Arm, Leg, Foot            |
| GL-703 | 3    | Medium Hand, Arm, Leg, Foot           |
| GL-704 | 4    | Large Hand, Arm, Leg, Foot            |
| GL-705 | 5    | Small Head, Shoulder, Thigh           |
| GL-706 | 5.5  | Medium Head, Shoulder, Thigh          |
| GL-707 | 6    | Large Head, Shoulder, Thigh           |
| GL-708 | 7    | Small Chest, Back, Perineum, Axilla   |
| GL-709 | 8    | Medium Chest, Back, Perineum, Axilla  |
| GL-710 | 9    | Large Chest, Back, Perineum, Axilla   |
| GL-711 | 10   | X-Large Chest, Back, Perineum, Axilla |
| GL-712 | 11   | Special Sizing                        |
| GL-713 | 12   | Special Sizing                        |
| GL-714 | 13   | Special Sizing                        |
| GL-715 | 14   | Special Sizing                        |

The exact size of net to be used is determined by the patient's weight, size and bulkiness of the dressing.

| N° cat. | Taille | Application                                   |
|---------|--------|-----------------------------------------------|
| GL-701  | 1      | Doigts, orteils, poignet                      |
| GL-702  | 2      | (Petit) main, bras, jambe, pied               |
| GL-703  | 3      | (Moyen) main, bras, jambe, pied               |
| GL-704  | 4      | (Grand) main, bras, jambe, pied               |
| GL-705  | 5      | (Petit) tête, épaule, cuisse                  |
| GL-706  | 5.5    | (Moyen) tête, épaule, cuisse                  |
| GL-707  | 6      | (Grand) tête, épaule, cuisse                  |
| GL-708  | 7      | (Petit) poitrine, dos, périnée, aisselle      |
| GL-709  | 8      | (Moyen) poitrine, dos, périnée, aisselle      |
| GL-710  | 9      | (Grand) poitrine, dos, périnée, aisselle      |
| GL-711  | 10     | (Très grand) poitrine, dos, périnée, aisselle |
| GL-712  | 11     | Grandeur spéciale                             |
| GL-713  | 12     | Grandeur spéciale                             |
| GL-714  | 13     | Grandeur spéciale                             |
| GL-715  | 14     | Grandeur spéciale                             |

La taille exacte du filet à utiliser est déterminée selon le poids du patient ainsi que la taille et la grosseur du pansement.